

## Aspire ALTITUDE Vertical Lift chair



- The ALTITUDE Vertical Lift chair is designed to accommodate varying user and facility requirements for enhanced comfort and efficient care provision.
- Integrated power source allows full chair functionality throughout the care environment for all-day independence.
- Intuitive control: Large backlit buttons make operation a simple process and enhance user independence.
- Comfort and pressure relief: Multi-stretch vapour permeable seat and backrest.
- The outer chair surface is upholstered in medical-grade vinyl for longevity and infection control.
- Assisted transfer functionality: Vertical lift and removable armrests facilitates safe sit-to-stand and lateral transfers, simplifying patient handling.
- Optional oncology armrests (CHP227710) wide forearm support enhance patient comfort and assists with effective administration of treatment such as blood transfusions.
- Optional arm rest table (CHP227720) allows for placement of personal items or medical instruments during treatment.

✓ Contracted DVA productZ

✓ Contract approved

### Support @ Home

**Function:** Sitting furniture

**Category:** Under advice

Click [here](#) for more information

### Compliance

ISO13934-1 / ISO13937-1 / BS6807 / BS7175

## Warranty Information

- 1 Year

Code	Display Name	Total Width	Total Height	Safe Working Load	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
CHP227700	Aspire ALTITUDE Vertical Lift	710 mm	1200 - 1340	150 kg	520 mm	510 mm	560 - 700 mm	560 - 700 mm	710 mm

Code	Display Name	Total Width	Total Height	Safe Working Load	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
	Chair		mm						
CHP227710	ALTITUDE Chair Oncology Armrests	-	-	-	-	-	-	-	-
CHP227720	ALTITUDE Chair Armrest Table	-	-	-	-	-	-	-	-